


Iron Deficiency

FOR SELF-TESTING USE



Feeling tired? Feeling run down?

The Iron Deficiency test is a simple at-home finger-prick blood test to quickly check your iron levels.

Iron is an important mineral your body needs to make healthy red blood cells. If you don't have enough iron, your body cannot make enough red blood cells, which could result in a condition called iron deficiency anaemia.

Knowing the reason for "feeling tired" is essential for your overall health: 12-15% of women and 2% of men are iron deficient. If iron deficiency is left unchecked, it can affect your brain, physical fitness and immune system.



Results in 5 Minutes



User friendly design



Easy-to-Read Results

Specification

Intended Use	Detects Iron levels at a cut-off concentration of 30ng/mL
Storage	2-30°C
Specimen Type	Whole Blood
Time to Result	5 minutes

Shelf Life	2 years
Sensitivity	91.3%
Accuracy	93.8%
Specificity	96.2%

Components

- 1 Test Cassette
- 1 Lancet
- 1 Capillary Dropper
- 1 Package Insert
- 1 Buffer
- 1 Alcohol Pad
- 1 Test Pack

Procedure



I. Pull off and dispose cap of lancet

II. Clean fingertip with alcohol pad

III. Press the lancet against the fingertip

IV. Massage the pricked fingertip to obtain blood drop



V. Put capillary dropper in contact with blood and draw the blood to fill line

VI. Put the collected blood into the sample well

VII. Add 1 drop of the buffer

VIII. Wait for the coloured line and read the results at 5 mins

Interpretation



NORMAL

Two distinct coloured lines appear Both T & C. Your iron levels are normal.



ABNORMAL

One line appears in C region. Your iron levels are low.



INVALID

The Control line fails to appear. Repeat the test.