



Title: **Aloclair PLUS Gel  
Leaflet**

Ref: **Sinclair/Aloclair/003**

Client: **Sinclair is Pharma**

Date: **06/07/12**

Size: **150mm x 210mm (Width/depth)**

Colours

Black  
(K)

Non-Printed Colours

Cutter/  
Fold

Giraffe Creative Technology Ltd., 35 Redsands Drive,  
Fulwood, Preston, Lancashire PR2 6GG (United Kingdom)  
Tel: +44 (0)1772 655900 - Fax: +44 (0)1772 700750  
email: info@giraffeld.com

# aloclair<sup>®</sup>

## PLUS

### ALOCLAIR<sup>®</sup> PLUS GEL - INFORMATION LEAFLET

Aloclair Plus Gel works by forming a protective film in the oral cavity that covers the exposed nerves endings of the lesion (aphthas) thus avoiding further irritation and reducing pain . The hyaluronic acid and the Aloe vera in the formulation support the natural healing of damaged tissues

Aloclair Plus Gel is indicated for aphthous ulcers, aphthous stomatitis, minor lesions, chafing and traumatic ulcers caused by braces and ill-fitting dentures.

#### Instructions for use:

Apply 1 or 2 drops of gel so as to cover the whole mouth ulcer or lesion. Avoid direct contact of the applicator with the lesion. Do not touch the lesion with the tongue for at least 2 minutes so as to allow for the protective film to form. Close the cap after each application. Use 3-4 times a day, or as needed and avoid eating or drinking for at least 1 hour after application .

#### Ingredients:

Aqua, polyvinylpyrrolidone (PVP), maltodextrin, propylene glycol, PEG-40 hydrogenated castor oil, xanthan gum, potassium sorbate, sodium benzoate, sodium hyaluronate, aroma, benzalkonium chloride, disodium EDTA , sodium saccharin, dipotassium glycyrrhizate, aloe barbadensis.

#### Cautions:

The use of Aloclair Plus Gel is contraindicated in patients with a known history of hypersensitivity to any of the ingredients.

#### Precautions and warnings:

Do not use the product if the package appears damaged .Store at room temperature. Do not freeze. Keep out of the reach of children.

If you frequently suffer from mouth ulcers seek medical advice.

**Contents:** 8ml

#### Distributed in the UK by:

Sinclair IS Pharma, Whitfield Court, 30-32 Whitfield Street, London W1T 2RQ.

#### Distributed in Ireland by:

Ocean Healthcare Ltd, Unit 6C, Nutgrove Office Park, Nutgrove Avenue, Rathfarnham, Dublin 14, Ireland.

#### Manufacturer:

Sinclair Pharma Srl, Viale Restelli 29, 20124 Milan, Italy.

Aloclair<sup>®</sup> is a trademark of Sinclair Pharmaceuticals Limited



0477

AB032901/1