

LYCLEAR®

Eradicates head lice & eggs

EFFECTIVE in

15 min

CLINICALLY TESTED



SPRAY
+COMB

Read the leaflet carefully before commencing the treatment.

About head lice

The head louse is a tiny wingless insect that lives exclusively in hair and feeds on blood. It moves around using its 3 pairs of feet that allow it to attach itself to hair.

An adult louse is 3-4 mm long.

The female lives for 2 to 3 months and produces 200 to 300 eggs in her lifetime.

They hatch after 7 days and become adults in 10 days.

When feeding, lice inject saliva into the scalp, which can cause a hypersensitivity reaction resulting in severe itching and irritation.

Diagnosis is confirmed by thorough examination of the hair, looking for lice, or eggs which are easier to see.

Precautions before use

Do not use this product on:

- people allergic to any of the ingredients
- children under 2 years old.

Pregnancy and breastfeeding

There is no counter indication for pregnant women to use Lyclear Spray.

However, it is recommended to always seek the advice of your doctor or pharmacist before using a product when pregnant or breast-feeding.

Warnings

For external use only. Consult a doctor in case of accidental ingestion.

Avoid all contact with the eyes and mucous membranes (mouth, nose). In case of contact with the eyes, rinse thoroughly with water.

Don't use in the case of allergy to one of the ingredients. In the event of a rash breaking out, an irritation, or hyper-sensitivity stop the treatment immediately, and wash the hair with a regular shampoo.

Possible side effects

Itching and irritation around the eyes. Flaking of the scalp. Contact irritation and redness of the scalp, ears or neck.

Storage

Store at room temperature

Do not use after expiry date displayed on bottle

Keep out of reach of children.

About Lyclear Spray

Lyclear Spray is clinically tested and proven to be an effective treatment for head lice.

Its formula of mineral oil and dimethicone has a double physical action based on suffocation and dehydration which gives the maximum result.

Because it contains no pesticides, resistance is not possible.

Composition

Mineral oil and dimethicone.

LYCLEAR®

How to use Lyclear?

All you need to do is follow the step-by-step instructions below.

Step No.1: check for the presence of head lice

Use the enclosed fine toothed Lyclear anti-lice comb to check thoroughly the hair of each family member, preferably in daylight.

Treat all those who have got head lice, the same day, in order to avoid re-infestation.

1) It is easier to check dry hair.

2) First of all, brush to remove tangles, then check the hair lock by lock.

3) Check the hair roots, behind the ears and on the nape of the neck, it is often in these areas that lice and eggs are found.

After each use, carefully clean the comb in lukewarm, soapy water, so as to avoid any further infestation.

Step No.2: apply Lyclear Spray

Hold the spray approximately 10 cm from the hair and carefully spray a sufficient quantity of Lyclear solution evenly onto dry hair. Be sure to cover all the hair and scalp.

Massage well from the roots to tips paying particular attention to the neck and the area behind the ears. Leave it on for 15 minutes.

Step No.3: wash the hair

Wash the hair with a normal shampoo, and be sure to rinse thoroughly. Wash again if necessary.

Step No.4: remove the lice & eggs

Before drying the hair, using the enclosed fine toothed Lyclear anti-lice comb to remove lice and eggs.

1) Comb wet hair with the fine toothed anti-lice comb.

2) Divide the hair into 4 sections: from the forehead to the nape of the neck then from one ear to the other. For long hair, you can braid each section.

3) Treat each section by isolating small locks of hair so that the head is visible.

4) Comb each lock of hair starting as close as possible to the head.

5) Using the comb, remove each egg that is found.

6) Each lock of hair that is egg free can be put to one side.

7) In order to avoid any re-infestation, clean the comb after each lock with a cotton wool pad or ball and always comb from the top of the head to the bottom.

Step No.5: repeat the treatment 7 days later

Repeat the treatment 7 days later repeating steps 2 to 4.

This will ensure that the life cycle of the lice is completely broken.

Any egg that could have hatched since the first treatment can also be eliminated.

There is no risk in repeating the treatment several times.

Step No.6: decontaminate the environment

Decontaminate clothes, bed linen and towels and other items (brushes, combs etc.) that were in direct contact with the head by washing them at 60°C. Carefully vacuum clean the house and car, then throw away the vacuum bag.

Important note:

Head lice cannot develop resistance to Lyclear Spray, as it has a non-chemical mode of action.

However it is possible that people who have been treated can get lice again from other untreated individuals (for example at school). In this event, you can use Lyclear Spray again at no risk, it will still be as effective.

Manufacturer:

Chefaro Ireland Ltd., Farnham Drive, Finglas, Dublin 11, Ireland.

Distributor:

Chefaro UK Ltd., Hamilton House, Mabledon Place, Bloomsbury, London, WC1H 9BB, United Kingdom.

Chefaro Ireland Ltd., Farnham Drive, Finglas, Dublin 11, Ireland.

www.headliceadvice.net (UK)

www.headliceadvice.ie (Ireland)

CE

twoscompany

Design Studio Ltd

Client: Chefaro UK

Project: Lyclear Spray Leaflet


Contact: Elva O'Connor

Job No: 5306

Size: 145 x 135mm

Date: 26/03/2010

Colours:



Address: Victorian House, 29 Ware Road, Hertford, Herts SG13 7EB
E: info@twoscompanydesign.com
T: 01992 586587